

NEW

Refrigerator

Blood bank or pharmaceutical refrigerator, 0-10°C, drawer- or wire shelves



BR Blood bank, drawer shelves PR Pharm. refrigerator, wire shelves

Ideal for:

- blood bank- & pharmaceutical refrigeration

Features:

- CFC free refrigerant (R-404A)
- intelligent control of refrigerator for preventing overload of compressor
- built-in battery for power failure alarm
- automatic condensate removal, no drain line required
- wire shelves for PR-models (6 or 12 pcs. acc. to model)
- drawer shelves for BR-model (6 or 12 pcs. acc. to model)
- stainless steel interior and powder-coated steel body
- optional built-in chart recorder & remote alarm system
- double fold tempered glass window
- castors for easy movement
- CE certified and unique serial number for tracing

Safety mechanism:

- overheat and over-current protection
- sensor error detection
- door locking system
- power failure alarm

Controller:

- *digital PID controller*



Model	Blood bank refrigerator		Pharmaceutical refrigerator	
	BR-600	BR-1000	PR-600	PR-1000
Capacity	620 l	1140 l	620 l	1140 l
Temperature range & sensor	0°C - +10°C (at 30°C room temperature & 30% humidity), PT 100 sensor			
Temperature resolution	control: 1.0°C, display: 1.0°C			
Compressor	250 W			
Evaporator	auto exchange dual evaporator			
Defrost & circulation	automatic defrost and melt water evaporation, dual turbo fan forced-air circulation			
Shelves	6 x drawer shelves	12x drawer shelves	5 x wire shelves	10x wire shelves
Internal dimensions (W x D x H)	620 x 700 x 1430 mm	1140 x 700 x 1430 mm	620 x 700 x 1430 mm	1140 x 700 x 1430 mm
External dimensions (W x D x H), net weight	740 x 805 x 1930 mm, 134 kg	1260 x 805 x 1930 mm, 208 kg	740 x 805 x 1930 mm, 134 kg	1260 x 805 x 1930 mm, 208 kg
Power supply	1 Phase AC 120V, 60 HZ or AC 230V, 50/60 Hz			
Power consumption	650 W	750 W	650 W	750 W
Order number 230V	DH.RefB0600	DH.RefB1000	DH.RefP0600	DH.RefP1000
Order number 120V	DH.RefB0601	DH.RefB1001	DH.RefP0601	DH.RefP1001

Accessories: Suitable shelves can be found on page 96