












PIPETTE TIPS

Cat. No.	Volume	Length	Type	Description	Packaging
					
BOE 1-300-81-0	0,1-10 µl	31,2 mm	A	Top Line Natural Tip	bag / 1000 pcs. / carton 20.000 pcs.
BOE 1-300-96-0	0,1-10 µl		A	Top Line Natural Tip	96-tip Rackbox / pack 10 racks
BOE 1-301-96-0	0,1-10 µl		A	Top Line Natural Tip, sterile	96-tip Rackbox / pack 10 racks
BOE 2-301-96-0	0,1-10 µl		A	Top Line Filter Tip, sterile	96-tip Rackbox / pack 10 racks
					
BOE 03.033.0100	0,1-20 µl	40,5 mm	B	Standard Tip, Natural, long type	bag / 1000 pcs. / carton / 10.000 pcs.
					
BOE 1-350-81-0	0,5-10/20 µl	45,7 mm	C	Top Line Natural Tip, long type	bag / 1000 pcs. / carton / 20.000 pcs.
BOE 1-350-96-0	0,5-10/20 µl		C	Top Line Natural Tip, long type	96-tip Rackbox / pack 10 racks
BOE 1-351-96-0	0,5-10/20 µl		C	Top Line Natural Tip, long type, sterile	96-tip Rackbox / pack 10 racks
					
BOE 03.013.0103	2- 200 µl	52,6 mm	D	Standard Tip, yellow	bag / 1000 pcs. / carton / 25.000 pcs.
					
BOE 1-410-81-1	1-200 µl	53,4 mm	E	Top Line Yellow Tip, bevelled, graduated	bag / 1000 pcs. / carton / 35.000 pcs.
BOE 1-410-96-1	1-200 µl		E	Top Line Yellow Tip, bevelled graduated, nonsterile	96-tip Rackbox / pack 10 racks
BOE 1-411-96-0	1-200 µl		E	Top Line Natural Tip, bevelled graduated, sterile	96-tip Rackbox / pack 10 racks
					
BOE 2-416-96-0	1-20/30 µl		E	Top Line Filter Tip, natural, sterile	96-tip Rackbox / pack 10 racks
BOE 2-412-96-0	1-100 µl		E	Top Line Filter Tip, natural, sterile	96-tip Rackbox / pack 10 racks
BOE 2-413-96-0	1-200 µl		E	Top Line Filter Tip, natural, sterile	96-tip Rackbox / pack 10 racks
BOE 2-415-96-0	5-300 µl		E	Top Line Filter Tip, natural, sterile	96-tip Rackbox / pack 10 racks
					
BOE 732031	5-300 µl	52 mm	F	Natural Tip	box / 1000 pcs. / carton / 10.000 pcs.
					
BOE 03.023.0101	50-1000 µl	70,5 mm	G	Standard Tip, blue	bag / 1000 pcs. / carton / 10.000 pcs.
					
BOE 1-209-06-2	100-1000 µl	70,1 mm	H	Top Line Tip, clear blue	60-tip Rackbox / pack 10 racks
BOE 1-207-06-2	100-1000 µl		H	Top Line Tip, clear blue, sterile	60-tip Rackbox / pack 10 racks
BOE 2-209-06-2	100-1000 µl		H	Filter Tip, clear blue, sterile	60-tip Rackbox / pack 10 racks
					
BOE 1-600-80-2	100-1000 µl	84,3 mm	I	Top Line Tip, blue, bevelled, graduated	bag / 500 pcs. / carton / 10.000 pcs.
BOE 1-601-96-2	100-1000 µl		I	Top Line Tip, blue, bevelled graduated	96-tip Rackbox / pack 10 racks
BOE 1-601-96-0	100-1000 µl		I	Top Line Tip, clear, bevelled graduated, sterile	96-tip Rackbox / pack 10 racks
					
BOE 2-601-96-0	100-1000 µl		I	Top Line Filter Tip, sterile	96-tip Rackbox / pack 10 racks

BOECO Top Line Tips are certified to be free of any detectable human DNA, DNase, RNase and PCR-inhibitors.

Cat.No.	Volume	Length	Type	Description	Packaging
---------	--------	--------	------	-------------	-----------



BOE 9600650	1 - 5 ml	126 mm	J	Natural Tip, for SA Series	carton / 300 pcs.
-------------	----------	--------	---	----------------------------	-------------------

This Tip is also suitable for our previous Standard and SP pipettes with 2000 μ l and 5000 μ l as well as BIOHIT[®] pipettes



BOE 702601	0,5 - 5 ml	160 mm	K	Natural Tip, for GP Series	bag / 200 pcs. / carton 1000 pcs.
------------	------------	--------	---	----------------------------	-----------------------------------

This Tip is also suitable for Thermo Fisher Scientific Finnpiquette[®] F1 and F2 of 5 ml



BOE 9600640	2 - 10 ml	150 mm	L	Natural Tip, for SA Series	bag / 100 pcs.
-------------	-----------	--------	---	----------------------------	----------------

This Tip is also suitable for Thermo Fisher Scientific Finnpiquette[®] F1 and F2 of 10 ml



BOE 702604	1 - 10 ml	156,5 mm	M	Natural Tip, for GP Series	bag / 100 pcs. / carton 1000 pcs.
------------	-----------	----------	---	----------------------------	-----------------------------------

This Tip is also suitable for the old SP Pipettes with 10 ml volume, as well as for Eppendorf[®] and Gilson[®] Pipettes of 10 ml



**NEW BOECO RACKBOXES
FOR 96 TIPS**